

# Dosis® pillbox reminds you 1 - 4 times daily when to take your medications.

### Clock-unit with a pillbox

The clock-unit is attached to a pillbox.



# Flexible compartments

in the container can be adjusted for the quantities of tablets. Four spacious compartments make it easy to put in and take out the pills.



reminds you to take your medication up to for times a day. Simple press of a button deactivates alarm when you have taken your pills. The pillbox is small (4 x 10 cm) and discreet with adjustable compartments.



# Large-digit LCD display

shows the time and indicates number of dosage times selected per day. If a dosage is missed, the display shows one or more crossed number and "MISSED PILL" twinkle to remind the patient and/or caregiver.



#### **Beeping Alarm**



### **Voice Alarm**



Visual Alarm



## Separate buttons

for setting hour - minute - alarms. The alarms can all be set at one time and the settings can be locked.